

Burton Medical UV Light | UV501



Stock Photo - Options May Vary

FACTORY NEW

The Burton UV501 Medical UV Light is a compact diagnostic device that delivers reliable UVA output at approximately 365 nm, enabling professionals to detect skin conditions, contaminants, and surface issues with confidence. Its durable construction and user-friendly design support fast, efficient examinations in medical, veterinary, and industrial environments.

The UV501 features a dual UV bulb system that provides steady illumination while maintaining energy efficiency. The lightweight handheld design allows easy handling during extended use and adapts to dermatology, forensic analysis, and inspection applications.

This dependable device combines durability, mobility, and precision to streamline workflow and enhance diagnostic accuracy. The simple on/off operation and long-lasting bulb lifespan ensure consistent, professional-grade performance across multiple clinical settings.

FEATURES

- Produces consistent UVA light for accurate examinations and inspections
- Supports dermatology, veterinary, and forensic applications
- Ensures long-lasting use with 6,000-hour bulb lifespan
- Delivers portable performance with lightweight handheld design
- Provides flexible installation with multiple mounting options
- Enhances workflow efficiency with simple on/off operation
- Offers durable construction for consistent, long-term reliability
- Adapts to multiple industries including medical and industrial inspection

SPECIFICATIONS

- Light Source: 2 UV fluorescent bulbs
- Wavelength: Approx. 365 nm (UVA)
- Bulb Life: Up to 6,000 hours
- Power Supply: 120V standard
- Design: Handheld, compact form factor
- Dimensions: Approx. 6.75" W x 3.75" H
- Mounting Options: Floor stand, edge clamp, or surface mount

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)