

Burton Medical UV Light with Magnifier | UV502



Stock Photo - Options May Vary

FACTORY NEW

The Burton Medical UV502 UV Light with Magnifier is a handheld diagnostic device that delivers consistent UVA output at approximately 365 nm, making it ideal for dermatology, veterinary, and clinical applications. Its ergonomic design ensures comfortable handling during extended use, enabling professionals to achieve clear vision and improved examination efficiency.

This UV exam light includes a high-quality magnifier with 2.25x magnification and a 5-diopter lens, allowing users to examine fine details with enhanced clarity. The dual UV bulbs offer long-term illumination with an average lifespan of up to 6000 hours, supporting durable, reliable performance.

The device combines advanced functionality with a user-friendly design featuring an on/off switch for quick operation. Its lightweight, portable structure supports convenient use in professional medical and clinical environments, delivering consistent results for skin evaluations and detailed inspections.

FEATURES

- Delivers consistent UVA output for accurate skin and surface examinations
- Enhances visibility with 2.25x magnification for detailed inspections
- Uses optically polished glass lens for superior clarity and durability
- Provides long-lasting performance with 6,000-hour bulb lifespan
- Supports dermatology, veterinary, and clinical diagnostic applications
- Ensures easy handling with lightweight, ergonomic handheld design
- Offers reliable illumination with dual UV bulb configuration
- Maintains consistent performance with stable 120V power input

SPECIFICATIONS

- Light Source: 2 UV fluorescent bulbs
- Wavelength: Approx. 365 nm (UVA)
- Magnification: 2.25x total power
- Lens Type: 5-diopter, optically ground glass
- Viewing Area: 8-inch diameter (20 cm)
- Bulb Life: Up to 6,000 hours
- Power Supply: 115-120V standard
- Design: Handheld configuration

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com