

# Burton Super Exam LED Ceiling Mount | SELEDC



*Stock Photo - Options May Vary*

## FACTORY NEW

The Burton SELEDC - Super Exam LED Ceiling Mount delivers powerful, focused illumination for medical exams and procedures. It delivers up to 35,000 lux at 1 m for clear visibility. Additionally, a CRI above 95 ensures accurate rendering of tissue and skin tones.

The 4000K color temperature provides clean, natural white light. As a result, clinicians achieve consistent lighting performance across general and specialty applications. Designed for precision, the Super Exam LED features a flexible, drift-free spring-balanced arm.

This design allows smooth positioning without unwanted movement. Moreover, the handle-mounted power and dimming controls improve workflow efficiency. Users can quickly adjust brightness from 100% to 0%.

## FEATURES

- High-performing energy-efficient LED module with 50,000-hour rated life
- Flexible drift-free spring-balanced arm
- Power and dimming switch positioned near the handle
- SteriTouch® coating inhibits the growth of bacteria
- 35,000 lux (3251 fc) at 39.37" (1 m)

## SPECIFICATIONS

- Illuminance: 35,000 lux (3251 fc) at 39.37" (1 m)
- Bulb Life (hours): 50000
- Color temperature: 4000 K
- CRI >951
- Light sources: Light Emitting Diode, minimum 50,000-hour rated life

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)