

Luvis E100 Portable LED Exam Light



Stock Photo - Options May Vary

FACTORY NEW

The Luvis E100 Portable LED Exam Light is designed for medical examinations in smaller clinical spaces. It delivers extended operational life and reduced energy consumption compared to traditional lighting, lowering maintenance costs and operational complexity. The E100 features three selectable color temperature options (3,500K, 4,000K, and 4,500K) and is available in ceiling, mobile, and wall-mounted configurations.

FEATURES

- Excellent Illumination
- Ergonomic & Friendly Design
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 3 Selectable Color Temperature Options
- Intuitive and Easy Way of Control
- Intensity: 60,000 lux
- Color Rendering Index: 95

SPECIFICATIONS

- LEDs: 8
- CRI: 95
- Color Temperature: 3,500K /4,000K /4,500K
- Head Size: 9.84" x 9.68"
- Focus Size: 7" (Circle)/11.8" (Circle)
- Intensity: Max. 60,000 lux
- Power Consumption: 20W
- LED Life Span: 50,000hr

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)