

Burton UV Hand-Held Light with Magnifier | UV503



Stock Photo - Options May Vary

FACTORY NEW

The Burton UV503 UV Hand-Held Light with Magnifier is a versatile examination and surgical lighting tool designed for healthcare professionals. It features special fluorescent bulbs that produce UVA light at approximately 365 nm wavelength, combined with optional magnification capabilities. The hand-held design provides flexibility for a wide variety of clinical applications.

The magnifier models include optically ground and polished glass five-diopter lenses with 2.25x total power magnification and an 8-inch focal length. The fluorescent bulbs deliver an average lifespan of 6000 hours, providing reliable long-term performance with minimal heat generation.

Available in 120V standard configuration with 230/240V options, the UV503 is assembled in the USA and backed by a 5-year limited warranty. Models are offered with or without magnification to suit specific clinical needs.

FEATURES

- Special fluorescent bulbs produce UVA at approximately 365 nm wavelength
- Magnifier and white bulbs
- 4-light magnifier models include optically ground and polished glass five diopter magnifier lenses with 2.25 x total power magnification and focal length of 8" (20 cm)
- 6000-hour average bulb life
- 120V model standard; 230 / 240V models available
- Available without Magnifier
- 5-year limited warranty
- Assembled in the USA

SPECIFICATIONS

- Fluorescent bulbs
- Bulbs produce UVA at approximately 365 nm wavelength
- Optically ground and polished glass five-diopter magnifier lens with 2.25x total
- Power magnification: 2.25x
- A focal length of 8"
- Bulb Life: 6000-hour average
- 120V model standard; 230 / 240V models available
- 5-year limited warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)