

Burton AIM HI LED Exam Light | Portable Mount



Stock Photo - Options May Vary

FACTORY NEW

The Burton AIM HI LED Exam Light delivers powerful illumination with up to 100,000 lux at 1 m for exceptional visibility. A 4300K color temperature and CRI above 95 ensure accurate tissue representation. High-intensity LED modules reduce shadows while creating a wide, adjustable light pattern for clarity during procedures.

The system enhances efficiency with user-friendly controls and durable construction. Three-step dimming options allow quick brightness adjustments, while the rotating center handle enables smooth focus changes without interrupting workflow. The sealed design supports easy cleaning in demanding clinical settings.

The drift-free arm system allows precise positioning and consistent stability. The light supports ceiling, wall, or mobile configurations to fit diverse spaces. With up to 50,000 hours of LED life and a five-year warranty, it reduces maintenance costs and provides long-term value for healthcare facilities.

FEATURES

- High-intensity LED modules reduce shadows and improve visibility
- Large, adjustable light pattern for different procedures
- 3-step dimming (100% / 75% / 50%) with push-button controls
- Focus control via sterilizable center handle
- Smooth, sealed design for easy cleaning
- Drift-free arm system for stable positioning
- Flexible mounting and positioning options
- Long-life, energy-efficient LED system

SPECIFICATIONS

- 100,000 lux @ 1 m (AIM HI 100)
- Color temperature: 4300 K
- CRI (color rendering index): >95
- Light source: LED (3 modules)
- LED life: ~50,000 hours
- Power consumption: ~45 W (AIM HI 100)
- Focusing: adjustable via center handle
- Warranty: 5-year limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)