

# Burton AIM HI LED Exam Light | Single Ceiling Mounted



Stock Photo - Options May Vary

## FACTORY NEW

The Burton AIM HI LED Exam Light Single Ceiling delivers powerful, precise illumination with up to 100,000 lux at 1 m for exceptional visibility. A 4300K color temperature and CRI above 95 ensure accurate tissue representation. High-intensity LED modules reduce shadows while creating a wide, adjustable light pattern.

The Burton AIM HI LED enhances efficiency with user-friendly controls and durable construction. Three-step dimming options allow quick brightness adjustments, while the rotating center handle enables smooth focus changes without interrupting workflow. The LED system offers up to 50,000 hours of life, reducing maintenance costs. Its sleek, sealed design supports easy cleaning for demanding clinical settings.

The Burton AIM HI LED features a drift-free arm system for precise positioning and consistent stability. The light supports ceiling, wall, or mobile configurations to fit diverse spaces. Manufactured in the USA, it includes a five-year warranty for long-term value and reliable performance.

## FEATURES

- High-intensity LED modules reduce shadows and improve visibility
- Large, adjustable light pattern for different procedures
- 3-step dimming (100% / 75% / 50%) with push-button controls
- Focus control via sterilizable center handle
- Smooth, sealed design for easy cleaning
- Drift-free arm system for stable positioning
- Flexible mounting and positioning options
- Long-life, energy-efficient LED system

## SPECIFICATIONS

- Illuminance: 100,000 lux @ 1 m (AIM HI 100)
- Color temperature: 4300 K
- CRI (color rendering index): >95
- Light source: LED (3 modules)
- LED life: ~50,000 hours
- Power consumption: ~45 W (AIM HI 100)
- Focusing: adjustable via center handle
- Warranty: 5-year limited

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)