

Burton Super Exam LED Floorstand | SELEDFL



Stock Photo - Options May Vary

FACTORY NEW

The Burton SELEDFL Super Exam LED Floorstand delivers powerful, focused illumination for exams and procedures with up to 35,000 lux at 1 m for clear visibility. The 4000K color temperature provides clean, natural white light, and a CRI above 95 ensures accurate rendering of tissue and skin tones for consistent lighting performance across exam rooms and specialty settings.

Designed for mobility and precision, the floorstand features a stable rolling base for easy room-to-room transport and a flexible, drift-free spring-balanced arm that allows smooth, accurate positioning without interrupting workflow. Handle-mounted power and dimming controls improve efficiency, with a fully dimmable range from 100% to 0% supporting a wide range of procedures.

Built for durability, the Super Exam LED includes a 50,000-hour rated LED lifespan and energy-efficient 20-watt operation that reduces maintenance and lowers long-term operating costs. The antimicrobial SteriTouch® coating helps inhibit bacterial growth, and the smooth, sealed design allows fast and effective cleaning.

FEATURES

- High-performing energy-efficient LED module with 50,000-hour rated life
- Flexible drift-free spring-balanced arm
- Power and dimming switch positioned near the handle
- SteriTouch® coating inhibits the growth of bacteria
- 35,000 lux (3251 fc) at 39.37" (1 m)

SPECIFICATIONS

- Illuminance: 35,000 lux (3251 fc) at 39.37" (1 m)
- Color temperature: 4000 K1
- CRI >951
- Light sources: Light Emitting Diode, minimum 50,000-hour rated life
- Power: 20 Watts

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)