

# Burton Super Exam LED Wall Mount | SELEDW



Stock Photo - Options May Vary

## FACTORY NEW

The Burton SELEDW Super Exam LED Wall Mount delivers powerful, focused illumination for exams and procedures with up to 35,000 lux at 1 m for excellent visibility. The 4000K color temperature provides clean, natural white light, and a CRI above 95 ensures accurate rendering of tissue and skin tones for consistent, reliable lighting across exam rooms and specialty applications.

Designed for precision and ease of use, this wall-mounted exam light features a flexible, drift-free spring-balanced arm that allows smooth, stable positioning exactly where light is needed. Handle-mounted controls enable quick power and dimming adjustments during procedures, with a fully dimmable range from 100% to 0% supporting various lighting needs.

Built for durability, the Super Exam LED includes a 50,000-hour rated LED lifespan and energy-efficient 20-watt operation that reduces maintenance and lowers operating costs over time. The antimicrobial SteriTouch® coating helps inhibit bacterial growth, and the smooth, sealed design allows easy cleaning with disinfectants for long-term performance and hygiene.

## FEATURES

- High-performing energy-efficient LED module with 50,000-hour rated life
- Flexible drift-free spring-balanced arm
- Power and dimming switch positioned near the handle
- SteriTouch® coating inhibits the growth of bacteria
- 35,000 lux (3251 fc) at 39.37" (1 m)

## SPECIFICATIONS

- Illuminance: 35,000 lux (3251 fc) at 39.37" (1 m)
- Color temperature: 4000 K1
- CRI >951
- Light sources: Light Emitting Diode, minimum 50,000-hour rated life
- Power: 20 Watts

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)