

# Burton XenaLux Surgical Headlight | XN10H



Stock Photo - Options May Vary

## FACTORY NEW

The Burton XenaLux XN10H Surgical Headlamp delivers powerful, adjustable illumination for surgical and clinical procedures. This fiber-optic headlamp produces a bright circular beam that enhances visibility during ENT, general, and specialty surgeries, while the lightweight design fits comfortably for extended use.

The XN10H connects to a high-intensity xenon light source that supports consistent light output with minimal heat. The adjustable spot range lets surgeons tailor the beam size quickly, and the optical module and linkage arm allow for precise positioning during procedures.

The system arrives complete with the headlamp, illuminator, and necessary cables for clinical setup. A cool venting design supports reliable performance during longer cases, and an hours-in-use meter aids maintenance planning and lamp replacement tracking.

## FEATURES

- High-intensity surgical illumination
- Bright circular beam for focused visibility
- Adjustable spot range for precision control
- Lightweight headlamp for extended comfort
- Cool venting system reduces heat at source
- Hours-in-use meter for maintenance insight
- Ideal for ENT, general, and specialty surgery
- Trusted surgical lighting solution

## SPECIFICATIONS

- Light Source: 300 W xenon lamp
- Illuminance: 120,000+ lux at working distance
- Beam Type: Circular fiber-optic beam
- Spot Range: Adjustable diameter
- Cooling: Quiet, cool venting system
- Warranty: Limited 3 years (headlamp & illuminator)
- Cable: Standard 3.5 bundle cables included
- Compatibility: Fits standard surgical fiber-optic light sources

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)