

Clinton 7' L Floor Mounted Parallel Bars | 3-4007



Stock Photo - Options May Vary

FACTORY NEW

Clinton 7' L Floor Mounted Parallel Bars: 3-4007. Clinton Physical Therapy Equipment offers one of the best values in its class. It is built tough for the rigors of everyday physical therapy use and affordably priced for nearly every practice and facility.

View our full line of products, from treatment tables to staircases, parallel bars and mat platforms, and much more. Clinton Floor Mounted Parallel Bars feature parallel bars and uprights designed to mount directly to the floor for stability. The Width is adjustable, 1-1/2" one-piece, stainless steel handrails, and heavy-duty, square, telescoping, steel uprights.

Model 3-4007 also features easy height changes with spring-loaded plungers and numbered height indicator strips.

FEATURES

- Parallel bars and uprights designed to mount directly to the floor
- Width adjustable, 1-1/2" one-piece, stainless steel handrails
- Heavy-duty, square, telescoping, steel uprights
- Easy height changes with spring-loaded plungers and numbered height indicator strips
- The second height knob removes side to side play

SPECIFICATIONS

- Length: 7'
- Width Adjustment: 6"/side
- Height Range: 26" - 39"
- Weight Capacity: 400 lbs under normal physical therapy use when properly installed
- Warranty: 5 Year Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)