

# Clinton Adjustable One-Sided Staircase | 4-6501-36



Stock Photo - Options May Vary

## FACTORY NEW

Clinton Adjustable One-Sided Staircase: 4-6501-36. Clinton Physical Therapy Equipment offers one of the best values in its class. Furthermore, It is built tough for the rigors of everyday physical therapy use and affordably priced for nearly every practice and facility.

View our full line of products, from treatment tables to staircases, parallel bars, and mat platforms, and much more. Clinton Adjustable One-Sided Staircase features 1.5" stainless steel handrails, easy height adjustments with spring plungers, and numbered indicator strips. A second knob ensures rigidity.

The hardwood plywood base has a UV satin finish and assembles with steel fasteners. Slate Gray PVC edge material adds durability.

## FEATURES

- 1.5 diameter stainless steel handrails with end bumpers
- Fast easy height adjustments with self-indexing spring action plungers and numbered height indicator strips (no holes to line up or count)
- The second knob removes any side-to-side play for a rigid upright
- Hardwood plywood base with durable, UV satin, topcoat finish
- Assembled with steel fasteners (No glued wood dowels)

## SPECIFICATIONS

- Model: 4-6501-36
- Platform: 24"L x 36"W x 24"
- Each Step 10"L x 36"W x 6"H
- Height: 49" - 61.75"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)