

Clinton Adjustable Small Staircase | 4-6500-30

FACTORY NEW

Clinton Adjustable Small Staircase: 4-6500-30. Clinton Physical Therapy Equipment offers one of the best values in its class.

Furthermore, it is built tough for the rigors of everyday physical therapy use and affordably priced for nearly every practice and facility. View our full line of products, from treatment tables to staircases, parallel bars and mat platforms, and much more.

Clinton Adjustable Small Staircase Features: 1.5 diameter stainless steel handrails with end bumpers, fast easy height adjustments with self-indexing spring action plungers, and numbered height indicator strips (no holes to line up or count), and the second knob removes any side-to-side play for a rigid upright.



Stock Photo - Options May Vary

FEATURES

- 1.5 diameter stainless steel handrails with end bumpers
- Fast easy height adjustments with self-indexing spring action plungers and numbered height indicator strips (no holes to line up or count)
- The second knob removes any side-to-side play for a rigid upright
- Hardwood plywood base with durable, UV satin, topcoat finish
- Assembled with steel fasteners (No glued wood dowels)

SPECIFICATIONS

- Model: 4-6500-30
- Platform 24"L x 30"W x 12"H
- 2 Steps 10" 30"L x 6"W x
- 3 Steps 10" 30"L x 4"W x
- Overall Size 54"L x 34"W x 36.5" -50"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)