

# Clinton Adventure Series – Safari Treatment Table | 7940



*Stock Photo – Options May Vary*

## FACTORY NEW

The Clinton Adventure Series Safari Treatment Table is designed for school-age pediatric patients, featuring engaging photographic graphics that make the exam experience more comfortable and educational. The table combines durability with easy maintenance through its high-pressure laminate surfaces with permanently embedded graphics that resist peeling. The 2" firm foam padding is covered in premium, stain-resistant vinyl upholstery in Gunmetal.

The tabletop features scratch and mar-resistant surfaces with chip-resistant PVC edges and seamless, rounded corners for safety. The table includes 4-door storage with adjustable shelves and soft-close, Euro-style hinges for convenient supply organization. A paper dispenser is included as standard.

Available with or without a backrest, this treatment table is built for the demands of a busy pediatric practice. The durable construction and stain-resistant materials ensure long-term performance and easy cleaning between patients.

## FEATURES

- A Clinton's original theme
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high-pressure laminate and will not peel off
- 2" of firm foam padding
- 4 door storage with two adjustable shelves
- Scratch and mar resistant surfaces
- Premium, stain-resistant, woven, knit-backed vinyl upholstery in Gunmetal
- Soft-close, Euro-style door hinges

## SPECIFICATIONS

- Model Number: 7940
- Exam Table w/Adjustable Backrest Dimensions: 68"L x 27"W x 33"H
- Exam Table Flat Top (model number 7940-1) Dimensions: 68"L x 27"W x 33"H

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)