

Clinton Bariatric Parallel Bars | 3-2106



Stock Photo - Options May Vary

FACTORY NEW

Clinton Bariatric Parallel Bars: 3-2106. Clinton Physical Therapy Equipment offers one of the best values in its class.

View our full line of products, from treatment tables to staircases, parallel bars and mat platforms, and much more. The Clinton Bariatric Parallel Bars feature a 1-1/2" thick hardwood plywood base with an ultra-durable UV-coated satin finish.

Safety treads are at both ends of the platform. The middle upright includes a through pin lock for extra support. The bars have an extra-wide 48" width for bariatric patients.

Additionally, the width-adjustable 1-1/2" one-piece stainless steel handrails offer easy height adjustment. The spring-loaded plungers include a numbered height strip for convenience.

FEATURES

- Parallel bars are mounted to a 1-1/2" thick hardwood plywood base with ultra-durable UV-coated satin finish
- Safety tread at both ends of the platform
- Middle upright with through pin lock for extra support
- Extra-wide 48" width for bariatric patients
- Width adjustable 1-1/2" one-piece, stainless steel handrails

SPECIFICATIONS

- Handrail Construction: 1 1/2" one-piece, stainless steel
- Uprights: Heavy-duty, square telescoping steel
- Length: 10'
- Width Range: 20" - 34"
- Height Range: 26" - 39"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)