


Clinton Complete – General Exam Room: 6



Stock Photo – Options May Vary

FACTORY NEW

Clinton Complete – General Exam Room: 6. Clinton Complete exam rooms provide designer solutions with coordinated colors, flooring, and Clinton equipment for private practices or clinics.

Customize the space with Clinton tables and cabinets to fit your needs. The room comes with a 5-year warranty.

Clinton Complete Qualifying Purchase (3 pc group) features are listed below: Exam Table 8870: Family Practice Table, Gray Laminate, and Burgundy Upholstery.

FEATURES

- Full cabinet model table side storage access from either side
- two deep storage drawers on the end
- Two-position length adjustable stirrups
- Padded, pull-out leg rest
- Top sections are easily replaceable in the field

SPECIFICATIONS

- Model 8870: 58"/72" extended L x 27"W x 31"H
- Model 8042-99: 66"L x 18"W x 35"H
- Desk Height: 30"H
- Model 8242-9: 42" L x 12"W x 24"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)