

Clinton Exam Room Chrome Series Stool | 2102



Stock Photo - Options May Vary

FACTORY NEW

Clinton Exam Room Chrome Series Stool: Model 2102. Clinton Chrome Series stool is a classic. With sturdy construction and dependable spin lift, these stools offer practical, comfortable seating at very affordable prices.

Clinton's Chrome Series Stool is an excellent choice for multi-purpose office or clinic use. Available in 21 colors including 7 colors to coordinate with Clinton's Laminate Cabinetry or Tables. Clinton Chrome Base Stool Model 2102 features a 1" diameter, all-welded S-style frame, and a round foot ring.

It has firm polyfoam padding, a 4" thick seat for added comfort, and 2" rubber wheel ball bearing casters for smooth movement.

FEATURES

- Round foot ring
- 1" diameter, S-style, all-welded frame
- Firm polyfoam padding
- 4" thick seat for extra comfort
- 2" rubber wheel, ball bearing casters

SPECIFICATIONS

- Model: 2102
- Seat Diameter: 14"
- Base Diameter: 19.5"
- Height Range: 18.5" - 25.5"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)