

# Champion 435 Gas Lift Stool with Back



Stock Photo - Options May Vary

## FACTORY NEW

The Champion 435 Gas Lift Stool with Back is a professional-grade seating solution that provides essential lumbar support for long medical shifts. It helps clinicians maintain a healthy posture during delicate procedures and patient examinations. You can transition between workstations effortlessly thanks to the smooth-rolling casters.

Therefore, your team remains comfortable and focused on patient care. Adjusting the Champion 435 stool height is simple and fast. The integrated hand lever allows for immediate seat elevation or lowering.

Consequently, every staff member can find their perfect ergonomic fit in seconds. The 360-degree swivel feature offers a full range of motion within tight spaces.

## FEATURES

- **Pneumatic Height Adjustment:** Lever-operated gas lift provides seamless height changes for various users
- **Adjustable Backrest:** Offers dedicated support to the lower back to reduce fatigue during procedures
- **360-Degree Swivel:** Enables total rotation for easy access to tools and patient charts
- **Medical-Grade Vinyl:** Resists fluids and stains while standing up to harsh chemical cleaners
- **Five-Leg Stability:** Prevents tipping and ensures a secure foundation on different floor surfaces

## SPECIFICATIONS

- **Weight Capacity:** 250 lbs
- **Seat Height Range:** 17" - 22.25" (standard gas lift)
- **Seat Diameter:** 16"
- **Base Diameter:** 22"
- **Backrest Style:** Adjustable lumbar support

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)