

# Clinton Folding Parallel Bars | 3-3317

## FACTORY NEW

Clinton Folding Parallel Bars: 3-3317. Clinton Physical Therapy Equipment offers one of the best values in its class.



Stock Photo - Options May Vary

It is built tough for the rigors of everyday physical therapy use and affordably priced for nearly every practice and facility. View our full line of products, from treatment tables to staircases, parallel bars, and mat platforms, and much more.

Clinton Folding Parallel Bars Features: Width adjustable floor folding parallel bars with double stiffening rails Width adjustable and 1-1/2" one-piece, stainless steel handrails. Easy height adjustment with spring-loaded plungers and numbered height strip, heavy-duty, square telescoping steel uprights.

## FEATURES

- Width adjustable floor folding parallel bars with double stiffening rails
- Width adjustable, 1-1/2" one-piece, stainless steel handrails
- Easy height adjustment with spring-loaded plungers and numbered height strip
- Heavy-duty, square telescoping steel uprights
- The second height knob removes side to side play

## SPECIFICATIONS

- Length: 7'
- Width Range: 16" - 24"
- Height Range: 28" - 41"
- Weight Capacity: 350 Load capacity under normal physical therapy use
- Warranty: 5 Year Limited Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)