

Clinton Heavy Duty Lift and Load Weight Shelf | 5-3078



Stock Photo - Options May Vary

FACTORY NEW

Clinton Heavy Duty Lift and Load Weight Shelf: 5-3078. Clinton Physical Therapy Equipment offers one of the best values in its class. It's built tough to withstand the rigors of everyday physical therapy use and priced affordably for nearly every practice and facility.

From a full line of treatment tables to staircases, mat platforms, parallel bars, and much more, Clinton offers a variety of physical therapy equipment designed and engineered to meet a wide range of needs—all with one source of convenience. Clinton Heavy Duty Lift and Load Weight Shelf Features: The construction uses all hardwood plywood.

Double-thick uprights add durability, and bumper tee-molding finishes all exposed edges.

FEATURES

- All hardwood, plywood construction
- Double thick uprights add durability
- All exposed edges are finished with bumper tee-molding
- 3 adjustable shelves with easy lock-into-position stops
- 1" thick shelves are covered with durable, black, ABS plastic

SPECIFICATIONS

- Model: 5-3078
- Width: 30"
- Depth: 20"
- Height: 78"
- Upper shelf adjustment range: 55" to 59"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)