

Clinton Pediatric Series, Zoo Bus, Blood Drawing Station | 67020

FACTORY NEW

The Zoo Bus Blood Drawing Station is a pediatric-focused phlebotomy chair designed for clinical efficiency and patient comfort. It features padded surfaces throughout, including a flip arm, armrests, seat, backrest, and padded top section for infants. The chair incorporates durable no-peel graphics embedded in high-pressure laminate on both sides and includes levelers under each wheel for stability.



Stock Photo - Options May Vary

The station provides extensive storage with two compartments on each side and an additional storage area at the end, supporting staff efficiency in space-conscious facilities. A paper dispenser is integrated into the infant section for convenient access. The chair supports a 200 lbs. load capacity under normal use.

The Zoo Bus accommodates both pediatric and infant patients with its specialized design. Optional coordinating step stool is available to enhance accessibility for younger patients.

FEATURES

- Padded flip arm and armrests
- Padded seat and padded backrest
- Padded top for infants
- No-peel graphics embedded into high-pressure custom laminate on both sides
- Levelers under each wheel
- Large storage capacity
- 200 lbs. load capacity under normal use
- 2 Storage compartments on the side and one on the end

SPECIFICATIONS

- Model number 67020
- Seat Dimensions: 17"w x 14"d
- Flat Top Section Dimensions: 42"w x 25"d
- Overall Length: 63"
- Overall Width: 30"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)