

Clinton Pediatric/Cool Pals Phlebotomy Cart | 67231



Stock Photo - Options May Vary

FACTORY NEW

Clinton Pediatric/Cool Pals Phlebotomy Cart: 67231. For more than 3 decades, Clinton Industries has been producing quality medical equipment designed and engineered for superior clinical performance. Clinton is one of the broadest and best sources for a wide range of quality, practical medical equipment.

Our lines of medical equipment and accessories are easy to use, coordinated, and cost-effective. Clinton Pediatric/Cool Pals Phlebotomy Cart Features: Clinton Kids™ Cool Pals graphics embedded in laminate will not tear or peel. Capital Blue, powder-coated base with a lockable door on soft close hinges.

Two removable plastic storage trays on easy slide tracks. ClintonClean,™ one-piece, solid plastic top with no seamed edges to absorb pathogens.

FEATURES

- Clinton Kids™ Cool Pals graphics embedded in laminate will not tear or peel
- Capital Blue, powder-coated base
- A lockable door on soft close hinges
- Two, 3.5" x 11" x 16.5" removable plastic storage trays on easy slide tracks
- ClintonClean,™ one-piece, solid plastic top with no seamed edges to absorb pathogens

SPECIFICATIONS

- Model Number: 67231
- Dimensions: 20"W x 20.75"D (top)
- Base Diameter: 24"
- Height Range: 33" - 43"
- Cabinet Dimensions: 13.25"W x 18"D 10.5"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)