

Clinton Physical Therapy Mirrors



Stock Photo - Options May Vary

FACTORY NEW

Clinton Physical Therapy Mirrors are designed for durability and affordability in clinical settings. Available in Tri-Panel, Single Panel, Mobile, or Wall Mount configurations, these mirrors feature 1/4" thick distortion-free glass with ANSI 297 safety backing for reliable performance during patient treatment and assessment.

The mirrors are constructed with hardwood and plywood frames finished with almond laminate backing. Mobile models roll smoothly on 3" swivel casters and include a cuff weight rack on the back for convenient storage and organization of therapy equipment.

View angles adjust easily to accommodate different treatment needs and patient positions. All mirrors come with a 5-year limited warranty, providing long-term value for physical therapy practices and rehabilitation facilities.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)