

Clinton Platform Mounted Parallel Bars | 3-2012



Stock Photo - Options May Vary

FACTORY NEW

Clinton Platform Mounted Parallel Bars: 3-2012. Clinton Physical Therapy Equipment offers one of the best values in its class. It is built tough for the rigors of everyday physical therapy use and affordably priced for nearly every practice and facility.

View our full line of products, from treatment tables to staircases, parallel bars and mat platforms, and much more. Clinton Platform Mounted Parallel Bars Features: Parallel bars mounted to 1-1/2" thick hardwood plywood base with ultra-durable, UV-cured satin finish. Safety tread at both ends of the platform.

Width is adjustable, 1-1/2" one-piece, stainless steel handrails with Heavy-duty, square, telescoping, steel uprights.

FEATURES

- Parallel bars mounted to 1-1/2" thick hardwood plywood base with ultra-durable, UV-cured satin finish
- Safety tread at both ends of the platform
- Width adjustable, 1-1/2" one-piece, stainless steel handrails
- Heavy-duty, square, telescoping, steel uprights
- Easy height changes with spring-loaded plungers and numbered height indicator strips

SPECIFICATIONS

- Model: 3-2012
- Length: 12'
- Width Range: 18" - 28"
- Height Range: 26" - 39"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)