

Clinton Platform Mounted 7' Parallel Bars | 3-2007

FACTORY NEW

Clinton Platform Mounted 7' Parallel Bars: 3-2007. Clinton Physical Therapy Equipment offers one of the best values in its class.



It's built tough for everyday physical therapy use and priced affordably for nearly every practice and facility. View our full line of products, from treatment tables to staircases, parallel bars and mat platforms and much more.

Stock Photo - Options May Vary

FEATURES

- Parallel bars mounted to 1-1/2" thick hardwood plywood base with ultra-durable
- UV-cured satin finish Safety tread at both ends of the platform
- Width adjustable
- 1-1/2" one-piece, stainless steel handrails
- Heavy-duty, square, telescoping, steel uprights

SPECIFICATIONS

- Length: 7'
- Width Range: 18" - 28"
- Height Range: 26" - 39"
- Warranty: 5 Year Limited Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)