

Clinton Staircase Storage Island | 4-8020



Stock Photo - Options May Vary

FACTORY NEW

Clinton Staircase Storage Island: 4-8020. Clinton Physical Therapy Equipment offers one of the best values in its class. Clinton PT equipment is built tough for the rigors of everyday physical therapy use and affordably priced for nearly every practice and facility.

View our full line of products, from treatment tables to staircases, parallel bars and mat platforms and much more. Clinton Staircase Storage Island Features: All hardwood plywood construction staircase, small corner staircase for limited space, and requires a minimum 9' ceiling.

Comes with safety tread on all steps and platform. All laminate band rack and dumbbell and cuff strips and holds a load capacity of 500 lbs.

FEATURES

- All hardwood plywood construction staircase
- Small corner staircase for limited space
- Requires a minimum 9' ceiling
- Safety tread on all steps and platform
- All laminate band rack and dumbbell and cuff strips

SPECIFICATIONS

- Model: 4-8020
- Platform 24"L x 30"W x 24"H
- Each Step 10"L x 30"W x 6"H
- Overall Size 63"L x 39"W x 59.5"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)