

Clinton Upholstered Mat Platform with Adjustable Backrest | 241-68



Stock Photo - Options May Vary

FACTORY NEW

Clinton Upholstered Mat Platform with Adjustable Backrest: 241-68. Clinton Physical Therapy Equipment offers one of the best values in its class. It is built tough for the rigors of everyday physical therapy use and affordably priced for nearly every practice and facility.

View our full line of products, from treatment tables to staircases, parallel bars and mat platforms and much more. Clinton 241-68 Features: The 241-68 has a 6' x 8' table which has a 2-piece, split backrest with each side 36" wide. 6-position, all-steel, adjustable backrest (28" long) and the 6-leg design with center support adds stability, along with triple-bolted, solid hardwood corner legs and plywood frame.

Frame attached with steel angles.

FEATURES

- Features 6-position, all-steel, adjustable backrest (28" long)
- 6' x 8' table has 2-piece, split backrest with each side 36"
- 6-leg design with center support adds stability
- Triple-bolted, corner legs
- Solid hardwood legs and plywood frame

SPECIFICATIONS

- Model: 241-68
- Length: 8'
- Width: 6'
- Height: 18"
- Backrest-model 241-68 only (each side): 28"L x 36"W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)