

Clinton Wall Mounted Folding Parallel Bars | 3-3307



Stock Photo - Options May Vary

FACTORY NEW

Clinton Wall Mounted Folding Parallel Bars: 3-3307. Clinton Physical Therapy Equipment offers one of the best values in its class. It is built tough for the rigors of everyday physical therapy use and affordably priced for nearly every practice and facility.

View our full line of products, from treatment tables to staircases, parallel bars and mat platforms and much more.

Clinton Wall Mounted Folding Parallel Bars Features:

Stainless steel handrails, Requires only 14" of floor space when stored and Folds up against the wall when not in use (folded height 92"). Rails are 1-1/2" one-piece, stainless steel handrails.

Easy height adjustment with spring-loaded plungers and numbered height strip on front legs.

FEATURES

- Stainless steel handrails
- Requires only 14" of floor space when stored
- Folds up against the wall when not in use (folded height 92")
- 1-1/2" one-piece, stainless steel handrails
- Easy height adjustment with spring-loaded plungers and numbered height strip on front legs

SPECIFICATIONS

- Length: 7'
- Width: 28"
- Height Range: 28" to 42"
- Folded Height: 92"
- Weight Capacity Under Normal Use: 400 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)