

Conmed Excalibur Plus PC – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Refurbished Conmed Excalibur Plus PC is a versatile electrocautery unit designed for general surgical applications. It delivers 300 watts of cutting power, 120 watts of coagulation, and 120 watts of bipolar capability. The unit features Aspen Return Pad Monitoring technology to help prevent patient site burns during procedures.

The system includes a digital display for easy operation and audible alarms to alert staff to important conditions. A monopolar footswitch is included as standard equipment. Optional bipolar footswitch and cart accessories are available for purchase to customize the unit to specific facility needs.

FEATURES

- Versatile, Powerful Unit
- 300W Cut/ 120W Coag/ 120 Bipolar
- Aspen Return Pad Monitoring
- Digital Display
- Audible Alarms
- Includes Monopolar Foot-switch

SPECIFICATIONS

- 300W Cut/ 120W Coag/ 120 Bipolar

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)