

Aspen Labs / Conmed Excalibur ESU – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Aspen Labs/Conmed Excalibur 1C is a versatile electrosurgical unit designed for general surgical applications. It delivers 300W cutting power, 120W coagulation, and 120W bipolar capability, making it suitable for a wide range of procedures. The unit features Aspen Return Pad Monitoring technology to help prevent patient site burns during use.

The Excalibur 1C includes a digital display for easy operation and audible alarms to alert users to important system conditions. It comes equipped with a monopolar foot switch for convenient hands-free control during procedures. The unit has been completely refurbished and tested to ensure reliable performance.

FEATURES

- 300W Cut/ 120W Coag/ 120 Bipolar
- Aspen Return Pad Monitoring Helps Prevent Site Burns
- Digital Display
- Audible Alarms
- Includes Monopolar Foot-switch
- Versatile, Powerful Unit

SPECIFICATIONS

- 300W Cut
- 120W Coag
- 120 Bipolar

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)