

Datascope Accutorr Plus Vital Signs Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Datascope Accutorr Plus is a multi-parameter vital signs monitor designed for screening in pre- and post-operative areas. The unit provides non-invasive blood pressure (NIBP), oxygen saturation (SpO2), and pulse rate monitoring with trending capabilities for enhanced clinical use.

This compact and portable monitor features a built-in pole clamp for flexible mounting and includes battery backup for continuous operation. The unit comes equipped with a new NIBP cuff, hose, new reusable SpO2 finger probe, and new battery. A temperature monitoring option is also available.

FEATURES

- Features NIBP, SpO2 (Datascope), and Pulse Rate
- Built-In Pole Mount Included
- Includes NIBP Cuff and Hose, SpO2 Finger Sensor
- Warranty: 6 Months - Parts and labor

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)