

Welch Allyn LXi Vital Signs Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Welch Allyn LXi Vital Signs Monitor – Refurbished. The refurbished Welch Allyn LXi Spot Vital Signs Monitor offers a reliable, convenient solution for fast, accurate patient monitoring.

Additionally, staff can easily move it from room to room, making it a good fit for busy clinics and hospitals. Furthermore, its compact size and light weight add mobility, while a built-in battery keeps it running without direct power.

Moreover, the LXi includes an easy-to-use carry handle, which simplifies transport between patients or locations. As a result, caregivers can deliver efficient, flexible care without losing accuracy or reliability.

FEATURES

- Features NIBP, SpO2 (Nellcor or Masimo), and Pulse Rate
- Optional Temp (SureTemp Oral/Rectal)
- 50 Reading Memory
- Lightweight Design
- Carry Handle

SPECIFICATIONS

- Dims: 9.7" H x 5.7" L x 8.75" W
- Weight: 4.5 lbs
- Battery: Rechargeable, Lead Acid
- Battery Voltage: 6V (2)
- Power Requirement: 120v

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)