

Welch Allyn LXi Vital Signs Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Welch Allyn LXi Vital Signs Monitor | Refurbished. The refurbished Welch Allyn LXi Spot Vital Signs Monitor offers a reliable and convenient solution for fast, accurate patient monitoring. It can easily be moved from room to room, making it ideal for busy clinics and hospital environments.

Moreover, its compact size and lightweight design enhance mobility, while the built-in battery ensures consistent performance even without direct power. Additionally, the LXi includes an easy-to-use carry handle, further simplifying transport between patients or locations.

As a result, caregivers can deliver efficient, flexible care without sacrificing accuracy or reliability.

FEATURES

- Features NIBP, SpO2 (Nellcor or Masimo), and Pulse Rate
- Optional Temp (SureTemp Oral/Rectal)
- 50 Reading Memory
- Lightweight Design
- Carry Handle

SPECIFICATIONS

- Dims: 9.7" H x 5.7" L x 8.75" W
- Weight: 4.5 lbs
- Battery: Rechargeable, Lead Acid
- Battery Voltage: 6V (2)
- Power Requirement: 120v

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)