

Edan iM3 Vital Signs Monitor – New



Stock Photo – Options May Vary

FACTORY NEW

The Edan iM3 Vital Signs Monitor is a portable vital signs monitoring device designed for fast-paced clinical environments. It features an 8-inch color touchscreen display, built-in WiFi connectivity, and support for hard key and knob operation. Standard parameters include SpO2, non-invasive blood pressure, and pulse rate.

The device offers multiple monitoring modes including continuous monitoring, spot checks, and a proprietary Round Mode optimized for rapid assessments. It includes a 360-degree alarm light, barcode scanner support, and connectivity options via USB, OTG, LAN, and DHCP. The lightweight, ultra-slim design makes it easy to transport between patient locations.

The iM3 integrates with Point Click Care through VitalsXChange and is compatible with Health o meter Professional Scales. Optional parameters are available including advanced SpO2, temperature measurement options, and alternative NIBP providers. The device runs on a rechargeable lithium-ion battery and operates on standard 120V power.

FEATURES

- SpO2
- Non-Invasive Blood Pressure
- Pulse Rate
- 8" Color TFT Touch Screen
- Built in WiFi

SPECIFICATIONS

- Dims: 6.3" W x 10.3" H x 6.5" W
- Weight: 5.5 lbs
- Battery: Rechargeable, Lithium-Ion
- Power Requirement: 120v

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)