

# Datascope Passport 2 Patient Monitor with CO2 - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Datascope Passport 2 Patient Monitor with CO2 is a refurbished solution offering excellent features and expandability. Built to meet the needs of both today and tomorrow, this monitor provides a cost-effective path for upgrades. As healthcare needs evolve, you can expand the capabilities of your Passport 2 system.

With future-ready features, this patient monitor provides flexibility while delivering reliable, high-quality patient care. The Datascope Passport 2 features a 10.4" high contrast color TFT display, offering exceptional visibility in any environment. It supports 3 to 6 waveforms, automatically expanding as needed for comprehensive monitoring.

The monitor adjusts large numerics automatically, ensuring that critical data is easily readable.

## FEATURES

- 10.4" high contrast color TFT display
- 3 to 6 waveforms with automatic expansion
- Auto-adjustable large numerics
- 3 or 5 lead ECG with ESU noise rejection
- Lead-Selectable Respiration

## SPECIFICATIONS

- Display: 10.4" high
- Temperature: 400 or 700 series YSI probes
- Graph and List Trending: Up to 120 patient time measurements
- OXY-CRG Trend: 6-minute for neonatal and pediatric applications
- Ports: Dual PCMCIA

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)