

Datex-Ohmeda Aestiva Anesthesia Machine - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Datex-Ohmeda Aestiva Anesthesia Machine - Refurbished. At the core of the Datex-Ohmeda Aestiva® 5 lies the SmartVent™ ventilator, known for its convenience, efficiency, and precision.

With both Pressure and Volume modes, the user-friendly SmartVent allows the Aestiva®/5 to effectively ventilate a wide range of patients, including neonates, trauma cases, compromised patients, and routine procedures. The Pressure Control mode lets you choose a target pressure and deliver the maximum tidal volume for the desired time and pressure selected.

The Volume Control mode delivers the tidal volume you set regardless of changes in the fresh gas flow and airway pressure - the pressure is selectable by the user pressure limit.

FEATURES

- SmartVent™ ventilator
- Integrated patient breathing system with circle module
- O2 and N2O gas delivery
- Pressure and Volume Controls enables the Aestiva®/5 to ventilate a broad range of patients: neonates, trauma, compromised and routine cases
- Compensation System: Optimizes low flow delivery techniques

SPECIFICATIONS

- Height: 53.4"
- Width: 36.6"
- Depth: 32.7"
- Approximate weight: 340 lbs.
- Work Surface Dimensions: 34.5"H x 18.5"W x 12.4"D

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)