

# GE Datex Aespire Anesthesia Machine - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The GE Datex Aespire Anesthesia Machine features the SmartVent ventilator, offering Pressure and Volume Control modes to ventilate a broad range of patients including neonates, trauma, compromised, and routine cases. The integrated breathing system simplifies cable management and improves patient safety by protecting components from damage and disconnection.

The SmartVent ventilator includes a compensation system that optimizes low-flow delivery techniques and automatically adjusts for changes in fresh gas flows, system leaks, lung compliance, and compression losses. This ensures more consistent delivery of set tidal volumes across varying clinical conditions.

The machine provides O<sub>2</sub> and N<sub>2</sub>O gas delivery with integrated O<sub>2</sub> and tidal volume monitoring. An open architecture allows integration with existing monitors and data management systems, while built-in service diagnostics reduce planned maintenance requirements and support costs.

## FEATURES

- SmartVent™ ventilator
- Integrated patient breathing system with circle module
- Pressure and Volume Controls enables the Aespire® to ventilate a broad range of patients: neonates, trauma, compromised and routine cases
- Compensation System: Optimizes low flow delivery techniques
- Smart sensors and switches: Provide uninterrupted communication between the breathing circuit and SmartVent™ ventilator

## SPECIFICATIONS

- Height: 53.4"
- Width: 36.6"
- Depth: 32.7"
- Approximate weight: 340 lbs.
- Warranty: 6 Months

*\*Brochure and images are a general overview of the product. Actual product options may vary.*