

GE Datex Aespire View Anesthesia Machine - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

GE Datex Aespire View Anesthesia Machine - Refurbished. The GE Aespire View offers advanced anesthesia in a compact system, ideal for busy operating rooms and surgical centers.

Additionally, it features the dependable GE SmartVent 7900 ventilator, offering precise control and ventilation performance. Furthermore, the unit also comes standard with GE's Advanced Breathing System (ABS), which supports low circuit volumes and low flow anesthesia, reducing agent consumption and operating costs.

Moreover, designed for reliability, ease of use, and safety, the Aespire View is a trusted choice for modern anesthesia care. GE collaborated with practitioners to design the Aespire View, focusing on features that truly impact daily anesthesia care and workflow.

FEATURES

- 7900 SmartVent ventilator
- 12" ventilation display
- Lightweight and Compact for easy maneuverability
- Electronic PEEP
- Smart alarms direct user to specific problems and affected parameters

SPECIFICATIONS

- Height: 53.5"
- Width: 29.5"
- Depth: 29.1"
- Approximate weight: 300 lbs.
- Warranty: 6 Months

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)