

Drager Narkomed 2B Anesthesia Machine - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The refurbished Drager Narkomed 2B Anesthesia Machine is durable and reliable. It includes O2 and N2O as standard, with an option for medical air.

The unit accommodates up to three vaporizers and features an AVE Ventilator with F-Style ascending bellows. Additionally, built-in airway pressure and oxygen monitoring ensure precise, reliable performance for various anesthesia needs.

Key features of the Narkomed 2B include clear, color-coded labeling for easy identification.

FEATURES

- 3 Gas (Air N2O, O2)
- 3 Vaporizers (19.1)
- AVE Electronic Ventilator
- Oxygen monitor and CO2 absorber and Scavenger Built-In, standard
- Oxygen Analyzer

SPECIFICATIONS

- Maximum height: 68"
- Maximum width: 40"
- Maximum depth: 25"
- Table-top height: 34"
- Table-top area: 271 square inches

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)