

# Drager Fabius GS Premium Anesthesia Machine - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Drager Fabius GS Premium Anesthesia Machine: Refurbished. The Drager GS Premium combines strong design with advanced ventilation features, making it ideal for most operating rooms. It mixes proven engineering with modern digital tech, offering great performance, value, and flexibility.

We fully refurbish and inspect this unit to ensure high reliability and quality. This Refurbished Fabius GS Premium anesthesia machine features a clear color screen and an easy Drager interface for simple operation.

It supports all key ventilation modes, making it flexible for different procedures. Additionally, the trolley is easy to move and has a central brake for stable positioning.

## FEATURES

- Affordable
- Reliable
- E-Vent piston ventilator
- Absorber
- PEEP

## SPECIFICATIONS

- Overall Dimensions: Approx. 30" x 52" x 33" without COSY, with 2 Vaporizer Mount
- Overall Dimensions: Approx. 40.7" x 52" x 33" with COSY and 3 Vaporizer Mount
- Weight: Approx. 296 lbs on base unit with COSY
- Operation Time with Fully Charged Batteries: Up to 120 minutes, minimum 45 minutes
- 6 months parts Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)