

Drager Apollo Anesthesia Machine – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Drager Apollo Anesthesia Machine: Refurbished. The refurbished Dräger Apollo Anesthesia Machine offers reliable performance, combining proven design with modern technology to meet the needs of today's operating rooms. It supports a wide range of patient needs and undergoes thorough inspection and refurbishment to ensure top performance and safety.

Moreover, the Apollo remains a trusted choice for medical facilities. The Apollo delivers ICU-quality ventilation with an electrically powered piston ventilator, eliminating the need for drive gas.

This simplifies setup and boosts efficiency. Its enhanced monitoring system gives real-time data, helping clinicians make informed decisions and ensuring reliability and precision during surgeries.

FEATURES

- Efficient and Ergonomic
- Reliable
- E-Vent plus
- Absorber
- PEEP

SPECIFICATIONS

- Overall Dimensions: Approx. 59" x 33.5" x 31.5"
- Weight: Approx. 365 lbs. (without vaporizers and gas cylinders)
- Power Backup: At least 30 minutes, typically 90 minutes, depending on ventilation procedures
- 6 months parts Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)