

Drive DeVilbiss 5 Liter Oxygen Concentrator | 525DS



Stock Photo - Options May Vary

FACTORY NEW

Drive DeVilbiss 5 Liter Oxygen Concentrator: 525DS. Since 2000, Drive Medical has been one of the fastest-growing global distributors of durable medical equipment.

The Company markets and distributes its products to customers throughout North America, Europe, Central America, South America, the Middle East, and Asia. Drive Medical carries a complete line of durable medical products, including respiratory equipment.

Drive DeVilbiss 5 Liter Oxygen Concentrator Features: Built upon a legacy of reliable, field-proven concentrator models, the DeVilbiss 5 Liter Oxygen Concentrator was designed with enhanced quality, durability, and simplicity in mind.

FEATURES

- Built upon a legacy of reliable, field-proven concentrator models
- Designed with enhanced quality, durability, and simplicity in mind
- Patented DeVilbiss Turn-Down Technology minimizes wear on internal components
- Turn-down technology reduces power consumption by 15% or more
- Turn-down technology also extends the life of the concentrator

SPECIFICATIONS

- Product Dimensions: 24.5"H x 13.5"W x 12"D
- Product Weight: 36 lbs.
- Power Consumption: 310-watt average, 275 watts @ 2.5 LPM and below
- Oxygen Percentage: 1-5 LPM, 93% +/- 3%
- Outlet Pressure: 8.5 psig

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)