

Ohio Care-E-Vac 3 Portable Aspirator Suction Unit – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Ohio Care-E-Vac 3 Portable Aspirator Suction Unit is built for fast response in compact, demanding clinical environments. Housed in a durable polycarbonate and ABS blend case, it's both tough and lightweight. Its small size and low weight make it perfect for code carts, transport, and mobile suction needs.

Furthermore, provides constant suction at 550 mmHg and 50 LPM — reliable power in critical moments. Moreover, designed for fast turnaround, it's easy to clean and simple to process between uses. The Care-E-Vac 3 is the smallest footprint in its class; it saves space without compromising performance.

Additionally, a permanently attached power cord keeps the unit ready to use — never lost or misplaced.

FEATURES

- Constant suction providing 550 mmHg and 50 LPM flow
- Easy to clean and process for quick turnarounds
- Compact and lightweight
- Smallest footprint in its class
- Permanently attached power cord cannot be lost or misplaced

SPECIFICATIONS

- Product Type: Aspirator
- Dimensions: 10.50" (L) x 8.5" (W) x 8.5" (H)
- Vacuum Range: 0-550 mmHg
- Free Air Flow: 50 L/min for 60 Minutes
- Pump: 12 VDC Oil-Less Diaphragm Type

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)