

SSCOR Duet 2314 Aspirator Suction Unit - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

SSCOR Duet 2314 Portable Aspirator Suction Unit - Refurbished is built for versatility and performance in critical care environments. Perfect for hospital crash carts, patient transport, surgical center backup, and long-term care facilities. This compact unit delivers strong, dependable suction where and when you need it most.

Additionally, operates on AC power and features a built-in battery backup for emergency use. Provides instant suction, even during power outages or patient transport scenarios. When fully charged, the battery offers approximately 45 minutes $\pm 10\%$ of run time.

The Duet 2314 is designed for quick setup and ease of use in urgent situations. Furthermore, each unit includes a canister, patient tube, power cord, and Quickstart Guide.

FEATURES

- AC suction unit
- Battery back up
- When fully charged at full capacity, 45 minutes $\pm 10\%$ running time
- >30 lpm clinical airflow
- Fully adjustable vacuum regulation between <50 and > 525 mmHg

SPECIFICATIONS

- Performance
- >30 lpm clinical airflow, <50 mmHg to > 525 mmHg vacuum
- Dimensions
- 9" H x 5.25 W x 17 L (22.9 cm H x 13.3cm W x 43.2cmL)
- Weight

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)