

Medco Aspirator IV Liposuction Aspirator



Stock Photo - Options May Vary

FACTORY NEW

The Medco Aspirator IV Liposuction Machine features dual positive displacement pumps that work together to relieve back pressure and draw vacuum directly from the patient. This setup ensures rapid, sustained deep vacuum performance during procedures. If one pump fails, the remaining pump continues working independently, ensuring continuous operation.

The aspirator incorporates the quietest positive displacement pumps available. Noise levels are minimized through 1-inch shock mounts that reduce vibration, a high-efficiency muffler, and an advanced insulation system. The Medco Aspirator IV delivers a comfortable environment for both patients and medical staff.

The unit is powered by standard AC electricity, while the foot pedal operates pneumatically, requiring no electrical power. This design fully complies with OSHA regulations, ensuring safety and reliability during use. The Medco Aspirator IV is made in the USA and meets the highest standards in medical equipment.

FEATURES

- High Efficiency: Dual-headed pumps ensure a deep vacuum is quickly reached and maintained.
- Quiet Operation: Equipped with a super quiet muffler and insulation system.
- Durability: Features oil-less piston vacuum motors for long-lasting performance.
- Versatility: Suitable for various medical applications beyond liposuction.
- Safety: Designed to prevent back pressure and ensure patient safety.
- Ease of Use: Includes a foot switch for hands-free operation.
- Reliability: If one pump fails, the other three continue to function independently.

SPECIFICATIONS

- Voltage: 110 or 220 V.
- Weight: 64 lb.
- Dimensions: 17" x 18" x 18".
- Max Vacuum: 29+ inches Hg.
- Max Flow: 7.5 CFM, which is approximately 212.38 LPM.
- Sound Level: 29 ± 3 db.
- Pump: 2 x 2 Piston.
- Horsepower: 2/3 hp.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)