


Medco Autofuse FR Super Flow Infiltration / Tumescant Pump



Stock Photo - Options May Vary

FACTORY NEW

The Medco Autofuse FR Super Flow Infiltration/Tumescant Pump delivers reliable infusion with forward and reverse flow capability for quick priming and purging of tubing. The intuitive design features a single-turn variable speed control knob for precise flow rate adjustments and a separate power switch that maintains speed settings when powering on or off, reducing procedural errors.

Dual foot pedal capability enables hands-free operation, allowing clinicians to focus on procedures without manual pump adjustment. This feature improves workflow efficiency and promotes consistent results during repetitive tasks.

Constructed from durable painted steel and ABS plastic, the pump is built to withstand demanding clinical environments. Its robust housing ensures long-term reliability and dependability for medical professionals.

FEATURES

- Forward and reverse flow motor allows easy prime/purge of the tubing in either direction.
- Simple operation with single-turn variable speed control.
- Separate power switch enables you to maintain the speed setting when turning the drive on/off.
- Dual foot pedal capability.
- Housing Material: Painted steel and ABS plastic.

SPECIFICATIONS

- Flow Range: 0.42 mL/min to 2,900 mL/min.
- Flow Rate Maximum: 2,900 mL/min.
- RPM Maximum: 600 RPM.
- Motor: 37 watts, 1/20 hp.
- Power: 0.8 amps, 50/60 Hz.
- Speed Control: \pm 5%.
- Weight: 16 lb.
- Dimensions: 9.1" (L) x 7.2" (W) x 5.3" (H) (231 mm x 184 mm x 134.5 mm).

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)