

# Gomco Optivac G180 Portable Aspirator Suction Unit - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Gomco Optivac G180 Portable Aspirator Suction Unit delivers high-performance suction in a compact, portable design. Ideal for general use, it operates on AC or DC power (110-230 V / 50-60 Hz) and is globally compatible. Additionally, its powerful motor provides over 30 lpm flow and a vacuum range of 25 to 550 mm Hg.

Furthermore, run the unit for up to three hours on a fully charged battery, even at maximum vacuum. Charge for two hours and get one hour of runtime; a full charge takes just eight hours. Use AC power at full performance, even if the battery is depleted.

Manufactured by Allied Healthcare, the Gomco Optivac G180 is engineered for reliability.

## FEATURES

- Run time operation on fully-charged battery for approximately three hours at maximum vacuum.
- 25 to 550 mm Hg vacuum range, coupled with greater than 30 lpm flow, provides powerful suction for faster aspiration.
- Will run at full performance on AC power even with a depleted battery.
- Will run for approximately one hour after charged for just two hours.
- Long-life rechargeable battery eliminates the need for frequent replacement.

## SPECIFICATIONS

- Vacuum Range: -25 to -550 mm Hg.
- Flow Rate: Greater than 30 lpm at open flow.
- Sound Level: Less than 60 dBA.
- Pump Type: Oil-less diaphragm.
- AC Operation: 110 to 230 volts, 50/60 Hz.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)