

# Drive VacuMax™ Go Suction Machine



Stock Photo - Options May Vary

## FACTORY NEW

Drive VacuMax™ Go Suction Machine: 18615. Easy-to-turn adjustable flow regulator knob conveniently placed directly below the gauge for easy use and control.

The model includes an 800cc canister, 6' tube, 10" tubing, hydrophobic filter, and elbow connector. It also comes with an internal rechargeable, replaceable battery that lasts up to 60 minutes per charge, and the manual.

Additional VacuMax™ Go Suction Machine Features: You can charge the 18615 while in use, and it requires less than 8 hours to reach a full charge. The Drive 18615 VacuMax™ Go is ideal for providing suction in the home or on the go.

## FEATURES

- Easy-to-turn adjustable flow regulator knob conveniently placed directly below gauge for easy use and control
- The model includes an 800 cc suction canister, 6' suction tube, 10" suction tubing, hydrophobic filter, elbow connector, rechargeable battery, carry bag, DC power cord for the car, and manual
- An internal rechargeable, replaceable battery that lasts up to 60 minutes per charge
- It may be charged while in use and requires less than 8 hours to reach a full charge
- Ideal for providing suction in the home or on the go

## SPECIFICATIONS

- Warranty on Suction Machine: Limited 2 Year
- Warranty on Battery: 6 months

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)