

Duraline WashChecks® Monitor Holder | C-WCH



Stock Photo - Options May Vary

FACTORY NEW

The Duraline WashChecks® Monitor Holder (C-WCH) is a reusable steel holder designed for secure placement and positioning of WashChecks® monitoring devices. This durable instrument is engineered for long-term use in sterilization environments and provides reliable support for process monitoring.

The holder accommodates both WashChecks® and WashChecks® U monitors, making it compatible with multiple monitor generations. Its hinged design allows for flexible positioning and measurement of steam impingement during autoclave cycles.

Constructed from steel, this holder is built to withstand repeated use and the demanding conditions of medical equipment sterilization. Sold individually, it serves as an essential accessory for proper monitor placement and operation.

FEATURES

- Holder is reusable
- Made to last forever (made of steel)
- Can be used with both the Wash-Checks® and Wash-Checks® U monitors
- The holder acts as a hinged instrument, measuring impingement

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)