

# DuraSonic .75 Gallon Digital Ultrasonic Cleaner | DS230D



Stock Photo - Options May Vary

## FACTORY NEW

DuraSonic DS230D .75 Gallon Digital Ultrasonic Cleaner. The new benchmark in bench-top cleaning! DuraSonic is the sound choice, now offering Ceramically Enhanced Transducers and a lifetime guarantee.

DuraSonic's Rugged Industrial-Strength Performance gives you consistent quality at a bench-top price. The driving force is the power source: while industrial stack transducers are known for reliability, bench-top cleaners use wafer-thin discs. DuraSonic now brings full-wave 45 kHz Ceramically Enhanced Industrial Transducers to the Benchtop line.

Each unit includes either 0-30 minute analog or 0-99 minute digital timers and thermostatically controlled heaters for efficiency. Compare features when evaluating bench-top cleaners—there is no comparison.

## FEATURES

- Digital model with timer, heat, power control, and degas
- Ceramically Enhanced Transducers
- All stainless steel enclosures resist damage from chemical spills and will not crack
- Sweep Frequency 45 kHz
- Adjustable, thermostatically controlled heaters (ambient to 80°C)

## SPECIFICATIONS

- Capacity: .75 Gallons
- Dimensions: 9.5" x 5.25" x 4"
- All Stainless Steel Enclosures
- Timer: 0 - 99 Minutes
- Lifetime Heater Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)